



**NICOTINE, COTININE, AND ANABASINE IN BLOOD,
PLASMA/SERUM, URINE, OR TISSUE BY LC-MS/MS OR GC-MS
CLEAN SCREEN[®] DAU EXTRACTION COLUMN**

Part #

ZSDAU020 – CLEAN SCREEN[®] DAU 200 mg, 10 mL Tube

SLPFPP100ID21-3UM – Selectra[®] PFPP HPLC Column, 100 x 2.1 mm, 3 μ m

1. PREPARE SAMPLE:

To 1 mL of 100 mM phosphate buffer (pH 6.0) add internal standards

Add 1 -2 mL of blood, plasma/ serum, urine, or 1 g (1:4) tissue homogenate

Mix/vortex and let stand for 5 minutes

Add 2 mL of 100 mM phosphate buffer (pH 6.0). Mix/vortex

Sample pH should be 6.0 ± 0.5 .

Adjust pH accordingly with 100 mM monobasic or dibasic sodium phosphate.

Centrifuge for 10 minutes at 2000 rpm and discard pellet

2. CONDITION CLEAN SCREEN[®] EXTRACTION COLUMN:

1 x 3 mL CH₃OH.

1 x 3 mL D.I. H₂O.

1 x 3 mL 100 mM phosphate buffer (pH 6.0).

NOTE: Aspirate at full vacuum or pressure

3. APPLY SAMPLE:

Load at 1 to 2 mL/minute.

4. WASH COLUMN:

1 x 3 mL D.I. H₂O.

1 x 2 mL 200 mM HCl

Dry column (5 minutes at full vacuum or pressure).

1 x 3 mL Methanol

Dry column (5 minutes at full vacuum or pressure).

5. ELUTE NICOTINE, COTININE, ANABASINE:

1 x 3 mL CH₂Cl₂/ IPA/ NH₄OH (78:20:2)

Collect eluate at 1 to 2 mL/minute.

NOTE: Prepare elution solvent daily.

Add IPA/ NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

6. DRY ELUATE:


Evaporate to dryness at < 40 °C.

7. RECONSTITUTE / DERIVATIZE

- **LC-MS/MS:** Reconstitute sample in 100 μ L of mobile phase
Inject 10 μ L.
- **GC-MS:** Dissolve residue in 100 μ L of Ethyl Acetate

INSTRUMENT CONDITIONS (LC-MS/MS):

CHROMATOGRAM SELECTRA® PFPP HPLC COLUMN



Analyte	MRM Transitions		Relative Retention Time (min)
	Q1	Q3	
1. Cotinine	117.2	80.1	2.42
Cotinine D ₃	180.2	101.2	-
2. Nicotine D ₄	167.2	136.1	3.03
3. Nicotine	163.2	132.2	3.06

PARAMETERS

Mobile Phase A: 0.1% Formic Acid in D.I. H₂O

Mobile Phase B: 0.1% Formic Acid in Methanol

Flow Rate: 0.3 mL/minute

Polarity: Positive

Reconstitute: 100 µL

Injection Volume: 10 µL

LC Column: Selectra® PFPP HPLC Column 100 x 2.1 mm 3 µm


Instrument: API 4000 Qtrap MS/MS with Agilent 1200 Binary Pump SL

Isocratic:

Time	%A	%B
0.00	80	20
10.00	STOP	

INSTRUMENT CONDITIONS (GC-MS):

CHROMATOGRAM



Analyte	Quantify Ion	Qualifier Ion 1	Qualifier Ion 2	Relative Retention Time (min)
1. Nicotine	84	133	162	5.71
Nicotine D ₄	88	137	166	-
2. Anabasine	84	105	133	6.32
3. Cotinine	98	119	176	7.12
Cotinine D ₄	101	122	179	-

PARAMETERS

GC/MS: Thermo ISQ Trace 1300

GC capillary column: 30 m x 0.25 mm (0.25 µm) TG-1MS

Injector: 1 µL Splitless, 250 °C

Oven temperature program: 50 °C (0.5) to 320 °C (30 °C/ minute): hold (5 minutes)

Carrier gas: Helium (1.2 mL/ minute)

MSD condition: Aux temperature: 280 °C, MS Source: 300 °C, MS Quad: 150 °C