



Method 523: Determination of Triazine Pesticides and their Degradates in Drinking Water by Gas Chromatography/Mass Spectrometry (GC/MS) Version 1.0

UCT Part Numbers: EC5232506 (250 mg GCB, 6 mL cartridge)

August 2012

Method Summary

This is a gas chromatography/mass spectrometry (GC/MS) method for the determination of triazine pesticides and their degradation products in finished drinking waters. Samples are pH adjusted, dechlorinated with ammonium acetate and protected from microbial degradation with 2-chloroacetamide during collection. Analytes are extracted from a **250 mL sample** using 250 mg carbon cartridges

The following compounds can be determined using this method:

Analyte	CASRN
Atrazine	1912-24-9
Atrazine-desethyl	6190-65-4
Atrazine-desethyl-desisopropyl	3397-62-4
Atrazine-desisopropyl	1007-28-9
Cyanazine	21725-46-2
Propazine	139-40-2
Simazine	122-34-9
Terbutylazine-desethyl	30125-63-4
Terbutylazine	5915-41-3
Prometon	1610-18-0
Prometryn	7287-19-6
Ametryn	834-12-8

1) Sample Preparation

- a. Allow samples to reach room temperature prior to extraction
- b. Add an aliquot of the Surrogate Primary Dilution Standards (PDS) to each sample
- c. Fortify Laboratory Fortified Blanks, Laboratory Fortified Sample Matrices, or LFSM Duplicates, with an appropriate volume of analyte PDS and the atrazine-desethyl-desisopropyl stock standard
- d. Cap and invert each sample several times to mix
- e. Proceed with sample extraction using SPE carbon cartridges

2) Cartridge Cleaning & Conditioning

- a. Set up extraction cartridges on the SPE vacuum manifold
- b. Using low vacuum (approximately 1 to 2 inches Hg), rinse each cartridge with two 6-mL aliquots of DCM drawing completely through
- c. Rinse each cartridge with a 6-mL aliquot of MeOH
- d. Draw MeOH to the top of the cartridge frit

Note: Do not let the cartridge dry after addition of MeOH

- e. Add a 6-mL aliquot of reagent water (RW) to the cartridge
- f. Draw RW to the top of the cartridge frit

3) Sample Extraction

- a. Add an additional 4 mL of RW to each cartridge
- b. Attach sample transfer lines to the cartridges. The additional volume prevents the SPE cartridge bed from going dry while the dead volume in the transfer lines is being filled
- c. Extract 250 mL of sample at a cartridge flow rate of 10 mL/minute
- d. Dry the cartridges under high vacuum for 10 seconds
- e. Release vacuum, then add a 0.25-mL aliquot of MeOH to each cartridge
- f. Draw the MeOH to waste, then dry cartridge under full vacuum for 10 minutes

4) Sample Elution

- a. Place 15-mL conical tubes into the manifold for collection
- b. Add 2 mL of EtOAc to the cartridge and elute dropwise
- c. Add 2 x 6-mL aliquots of 9:1 DCM/MeOH to cartridge
- d. Allow the cartridge beds to briefly soak in solvent, then draw the solvent through the cartridges
- e. Dry the eluate by passing it through approximately 3 grams of anhydrous Na₂SO₄ collecting it in a 40-mL centrifuge tube. Pre-rinse the Na₂SO₄ with a 1-mL aliquot of 3:1 DCM/EtOAc
- f. Rinse with 1-mL aliquot of 3:1 DCM/EtOAc collecting it in the centrifuge tube
- g. The dried extracts may be stored overnight in the 40-mL tubes at -10 °C
- h. Warm the 40-mL tubes to 35 °C in a water bath under a stream of N₂ and evaporate solvent to less than 1 mL but no less than 0.5 mL
- i. Transfer the concentrated eluate to 1-mL volumetric tubes
- j. Rinse the conical tube with a small volume of EtOAc, and transfer the rinse to the volumetric
- k. Add IS solution and adjust to volume
- l. Transfer the extracts to autosampler vials for analysis or store in a freezer ≤ -10 °C

Complete details at Office of Water (MLK 140) EPA Document No. 815-R-11-002 February 2011
<http://www.epa.gov/safewater/>

DCN-217280-240