



## COCAINE AND BENZOYLECGONINE IN URINE BY LC-MS/MS OR GC-MS STYRE SCREEN<sup>®</sup> DBX EXTRACTION COLUMN

Part #

SSDBX033 – STYRE SCREEN<sup>®</sup> DBX 30 mg, 3 mL Tube

SBSTFA-1-1 – SELECTRA-SIL<sup>®</sup> BSTFA w/ 1% TMCS

SLDA50ID21-5UM – SELECTRA<sup>®</sup> DA HPLC Column 50 x 2.1 mm, 5 $\mu$ m

### 1. PREPARE SAMPLE:

To 1 mL of urine add internal standard(s) and 300  $\mu$ L 100mM HCl.  
Mix/Vortex.

### 2. APPLY SAMPLE:

Load at 1 to 2 mL/minute.

### 3. WASH COLUMN:

1 x 1 mL D.I. H<sub>2</sub>O.

1 x 1 mL 100 mM HCl

1 x 1 mL CH<sub>3</sub>OH.

Dry column (10 minutes at full vacuum or pressure).

### 4. ELUTE COCAINE/METABOLITE:

2 x 0.5 mL CH<sub>2</sub>Cl<sub>2</sub>/ IPA /NH<sub>4</sub>OH (78:20:2)

Collect eluate at 1 to 2 mL/minute.

**NOTE:** Prepare elution solvent daily.

Add IPA /NH<sub>4</sub>OH, mix, then add CH<sub>2</sub>Cl<sub>2</sub> (pH 11-12).

### 5. DRY ELUATE:


Evaporate to dryness at < 40 °C.

### 6. RECONSTITUTE / DERIVATIZE:

- **LC-MS/MS:** Reconstitute sample in 100  $\mu$ L of mobile phase  
Inject 10  $\mu$ L.
- **GC-MS:** Dissolve residue in 50  $\mu$ L of Ethyl Acetate and 50  $\mu$ L BSTFA w/ 1%  
TMCS  
Overlay with N<sub>2</sub> and cap. Mix/vortex  
React 30 minutes at 70°C; Cool and inject 1  $\mu$ L

# INSTRUMENT CONDITIONS (LC-MS/MS):

## CHROMATOGRAM



Analyte	MRM Transitions		Relative Retention Time (min)
	Q1	Q3	
1. Benzoylcegonine D <sub>8</sub>	298.1	171.1	1.58
2. Benzoylcegonine	290.1	168.1	1.60
3. Cocaine	304.1	182.1	2.10

## PARAMETERS

**Mobile Phase A:** 0.1% Formic Acid in D.I. H<sub>2</sub>O

**Mobile Phase B:** 0.1% Formic Acid in Methanol

**Flow Rate:** 0.7 mL/minute

**Polarity:** Positive

**Reconstitute:** 100 µL

**Injection Volume:** 10 µL

**LC Column:** Selectra<sup>®</sup> DA HPLC Column 50 x 2.1 mm 5 µm


**Instrument:** API 3200 Qtrap MS/MS with Shimadzu Prominence UFLC

## Gradient:

Time	%A	%B
0.00	75	25
3.00	50	50
3.01	10	90
4.00	75	25
5.50	STOP	

## INSTRUMENT CONDITIONS (GC-MS):

### CHROMATOGRAM



Analyte	Quantify Ion	Qualifier Ion 1	Qualifier Ion 2	Relative Retention Time (min)
6. Cocaine	182	198	303	9.72
Cocaine D <sub>3</sub>	185	201	306	-
7. Benzoyecgonine TMS	240	256	361	9.94
Benzoyecgonine TMS D <sub>3</sub>	243	259	369	-

### PARAMETERS

**GC/MS:** Thermo ISQ Trace 1300

**GC capillary column:** 30 m x 0.25 mm (0.25 µm) TG-1MS

**Injector:** 1 µL Splitless, 250 °C

**Oven temperature program:** 70 °C (0.5) to 320 °C (25 °C/ minute): hold (2 minutes)

**Carrier gas:** Helium (1.2 mL/ minute)

**MSD condition:** Aux temperature: 280 °C, MS Source: 350 °C, MS Quad: 150 °C