

Development of Optimized Multiclass Clean-up Methods for LC-MS/MS Analysis of Mycotoxins in Multiple Food/Feed Matrices Incorporating a Novel SPE Column



Biotage®

Adam Senior¹, Geoff Davies¹, Claire Desbrow¹, Alan Edgington¹, Rhys Jones¹, Steve Jordan¹, Mats Leeman², Helen Lodder¹, Kerry Stephens¹, Lee Williams¹.

¹Biotage GB Limited, Distribution Way, Dyffryn Business Park, Ystrad Mynach, Cardiff, CF82 7TS, UK.
²MIP Technologies AB (a subsidiary of Biotage AB), Box 737, SE-220 07, Lund, Sweden.

Introduction

Mycotoxins are a group of toxic and often carcinogenic or genotoxic metabolites produced by several strains of fungi found on food crops worldwide. Mycotoxins have great potential to cause harm to humans, crops and farmed animals, example structures shown in **Figure 1**.




Figure 1. Typical Mycotoxin Structures – aflatoxin B₁, zearalenone, deoxynivalenol, ochratoxin A (clockwise from top left)

As a result, a wide range of food and feedstuff substrates require testing for mycotoxin contamination. The diversity of analyte structure and food substrate generates a significant analytical challenge. Traditionally, mycotoxins have been analyzed using multiple methods, each optimized for a single mycotoxin or group of closely related toxins. Multi-analyte approaches are being made increasingly possible by the adoption of liquid chromatography-tandem mass spectrometry (LC-MS/MS) based analyses. Due to the selective nature of LC-MS/MS, highly selective sample preparation techniques are often no longer essential to meet the residue limits required. However, appropriate sample preparation remains necessary in order to minimize matrix effects and maximize assay robustness.

Experimental

Reagents

Standards were purchased from LGC (Teddington, UK). Acetic acid (HAC), formic acid, ammonium hydroxide, ammonium acetate and toluene were purchased from Sigma-Aldrich Chemical Co. (Poole, UK). Chromatographic solvents were LC/MS grade and purchased from LGC/Promochem (Teddington, UK). 18.2 MΩ-cm water was drawn fresh daily from a Direct-Q 5 (Merck Millipore, Watford, UK)

Sample Preparation

Sample extraction: Optimized extraction conditions are tabulated below in **Table 1**. Briefly, a sub-sample of substrate was mixed with a specific volume of extraction solvent. The mixture was then extracted using the method and duration shown in **Table 1**.
Extraction work-up: The extract was centrifuged at 4000 x g for 10 min. Supernatant was transferred to a 50 mL tube and diluted with water to give a 10% final concentration of organic solvent. The diluted extract was centrifuged at 4000 x g for 10 min. Pre-treatment for small polar mycotoxins such as deoxynivalenol required separate pre-treatment 100% aqueous conditions.

Table 1. Substrate-Optimized Extraction Conditions

Substrate	Mass / g	Extraction Solvent	Solvent Vol. / mL	Extraction Method	Duration / min
Grain	5	50% acetonitrile (aq)	20	shake	30
Feed	5	1% HCOOH / 75% acetone (aq)	40	shake	30
Nuts	5	80% acetonitrile (aq)	20	shake	30
Chili (dried)	5	80% acetonitrile (aq)	20	Shake	30
Infant Formula (dried)	as mfr instructions	1% HCOOH (aq)	as mfr instructions	ultra-sonicate	20
Apple (juice)	1:1	1.0mM ammonium acetate pH 5 + adjusted to pH 5	20 (total)	mix	N/A

Solid Phase Extraction Clean-up: ISOLUTE® Myco SPE Columns, 60 mg / 3 mL (tabless); 150-0006-BG. **Table 2** below illustrates substrate-specific optimized SPE procedures.

Table 2. Substrate-Optimized ISOLUTE Myco SPE Procedures

Step	Grain & Feed	Nuts	Chili	Infant Cereal	Infant Formula	Juice (apple)
Condition (2 mL)	MeCN	MeCN	MeCN	MeCN	MeCN	MeCN
Equilibrate (2 mL)	H ₂ O	10mM NH ₄ OAc	H ₂ O	H ₂ O	H ₂ O	10mM NH ₄ OAc pH5
Load (3 mL)				(5 mL)	(5 mL)	(1 mL)
Wash 1 (2 mL)	H ₂ O	10mM NH ₄ OAc	H ₂ O	H ₂ O (5 mL)	H ₂ O (5 mL)	H ₂ O (3x 1mL) dry 5 min
Wash 2 (2 mL)	10% MeCN	10mM NH ₄ OAc / 10% MeCN	10% MeCN	10% MeCN (5 mL)	10% MeCN (5 mL)	Toluene (1 mL) dry 5 min
Elute 1 / Wash 3 (2 mL)	0.1% HCOOH / MeCN	0.1% HCOOH / MeCN	0.1% HCOOH / 40% MeCN	0.1% HCOOH / MeCN (5 mL)	0.1% HCOOH / MeCN (5 mL)	0.1% HCOOH / MeCN (1 mL)
Elute 2 (2 mL)	0.1% HCOOH / MeOH	0.1% HCOOH / MeOH	1% NH ₄ OH (conc) / MeOH	0.1% HCOOH / MeOH	0.1% HCOOH / MeCN	N/A

Clean-up for small polar mycotoxins such as deoxynivalenol required 100% aqueous load/wash conditions (details supplied on request).
SPE Work-up: The eluates were combined where necessary and dried under a stream of air at 45 °C for a period dependent on eluent composition. The eluates were reconstituted in 1 mL 0.1 % acetic acid in 1:1:8 MeOH:MeCN:H₂O (v/v/v) and syringe-filtered using a 0.2 µm PTFE membrane.

HPLC Conditions

Instrument: Waters Alliance 2795 (Waters Assoc., Milford MA, USA) or Shimadzu Nexera (Shimadzu Europa GmbH, Duisburg, Germany).
Patulin HPLC Column: Kinetex Phenyl-Hexyl 50 x 2.0 mm id, 2.6 µm (Phenomenex, Cheshire UK).
Mycotoxins HPLC Column: Kinetex C18-XB 50 x 2.0 mm id, 2.6 µm (Phenomenex, Cheshire UK).

Mass Spectrometry

Instrument: Ultima Pt triple quadrupole mass spectrometer (Waters Assoc., Manchester, UK) with electrospray interface or Triple Quad 5500 mass spectrometer (AB Sciex, Framingham, US) with a Turbo-V electrospray interface.

ESI Analysis: Positive and/or negative pseudomolecular ions were fragmented to their most abundant product ions for MRM (Waters) or sMRM (AB Sciex) experiments.

Chromatography and mass spectrometry experimental details can be supplied on request.

Results

Resin Development

Sample pre-treatment conditions were established during resin development. We found that 50% MeCN gave good recoveries for most classes of mycotoxins. However, for mycotoxins with significant polar character (e.g. deoxynivalenol): extraction, SPE equilibration and load steps required 100% aqueous conditions. SPE conditions were established during resin development. A 60 mg resin bed binds 100% of a suite of mycotoxins when loading up to 10 mL of 200 µg kg⁻¹ spiked flour extract (data not shown). The optimum % organic in the wash solvent was determined to be 10% as demonstrated in **Figure 2**. This can




Figure 2. Cumulative elution profile of mycotoxins varying wash procedure

be increased to 20% where recovery of aflatoxins is not required. Equilibration/wash conditions were investigated. Optimum SPE performance was obtained using weak solvents, typically H₂O. Low concentration ammonium acetate was found to improve recoveries where the substrate contained lipids or had increased ionic strength.

Reconstitution

Recoveries were stabilized using 0.5% HAC in preference to 0.1% (and HCOOH) as demonstrated in **Figure 3**; 1:1:8 MeOH:MeCN:H₂O (v/v/v) gave the best results, 20% MeOH was an acceptable alternative.

Recoveries

ISOLUTE Myco Analyte Recovery at EU MRL




Figure 4. Recovery of a suite of spiked mycotoxins from grains using ISOLUTE Myco (2 µg kg⁻¹ aflatoxins, 20 µg kg⁻¹ others) between 70% and 110% of spiked.

The recovery of mycotoxins from ground soya or animal feed using the optimized method from **Table 2** is demonstrated in **Figure 5**. 19/20 analyte recoveries were between 85% and 109%.

Recovery of Mycotoxins from Peanut and Soybean Nut Paste using an Optimized Method at Low Spike Levels



Figure 6. Recovery of mycotoxins from two nut matrices using a SPE method optimized for ISOLUTE Myco (1.6 µg kg⁻¹ spike)

non-polar organic solvent washes for interference removal gave no observable benefit to the method. The recovery of mycotoxins using the optimized method is demonstrated in **Figure 6** and was 89% to 105% (113% for AG1 in peanuts) at a spike of 1.6 µg kg⁻¹.




Figure 3. Comparison of reconstitution conditions on mycotoxin suite recoveries from ground wheat

Figures 4 through 10

demonstrate SPE recoveries at EU MRLs using the optimized methods shown in **Table 2**. The recovery of a suite of mycotoxins from differing grain substrates is demonstrated in **Figure 4**. With the exception of DON, recoveries were




Figure 8. Recovery of mycotoxins from infant cereal using a SPE method optimized for ISOLUTE Myco (0.2 to 2 µg kg⁻¹ spike)

Where substrates had high lipid content (e.g. nuts), recoveries of ochratoxin A were optimized by increasing MeCN concentration to 80% and changing the equilibration/wash solvent to 10 mM ammonium acetate. The use of hexane and salt as defatting agents; and the use of non-polar organic solvent washes for interference removal gave no observable benefit to the method. The recovery of mycotoxins using the optimized method is demonstrated in **Figure 6** and was 89% to 105% (113% for AG1 in peanuts) at a spike of 1.6 µg kg⁻¹.

Pigment and phenolic co-extractives in ground chili gave rise to low S/N and high suppression.

The effects were mitigated by reducing the % organic in elute 1 and adding 1% NH₄OH to elute 2. Recoveries are demonstrated in **Figure 7**.

Figure 7. Recovery of mycotoxins from ground chili using a SPE method optimized for ISOLUTE Myco (2 to 10 µg kg⁻¹ spike)

Infant cereal presented the additional challenge of very low pbp MRLs. We were able to meet these for the most part as demonstrated in **Figure 8**. However, aflatoxin B₁ LOQ was above the EU MRL and aflatoxin G₂ recovery was below the EU criterion due to matrix effects observed at the same retention time.




Figure 5. Low-level recovery of mycotoxins spiked into animal feed or ground soya using the optimized ISOLUTE Myco Method (5 µg kg⁻¹ aflatoxins, 100 µg kg⁻¹ Z/F/D/T spike)

Determination of AM1 at an adequate level as demonstrated in **Figure 9** required the use of AB₂ as a surrogate internal standard. The use of isotopically labelled AM1 is expected to give a more robust method.

The benefit of alternative wash strategies for widely differing substrates is demonstrated in **Figure 10**. Patulin recovery is lower using a toluene interference wash. However, signal is increased (not shown) and suppression is reduced to tolerable levels. Reduced signal suppression brings a concomitant improvement to method LOQ. Validation data for the presented recoveries above are available on request.




Figure 10. Recovery and signal suppression of patulin with and without a toluene interference wash (100 ng mL⁻¹ patulin spike)

Conclusions

- The optimized selectivity of the ISOLUTE Myco polymer allows for simultaneous extraction/measurement of multiple mycotoxins.
- The novel polymer used in ISOLUTE Myco is able to recover between 70% and 110% of spike for a majority of mycotoxins from a wide range of substrates using simple catch-and-release methods to below most EU MRLs.
- ISOLUTE Myco demonstrates sufficient sample clean-up to facilitate robust determination of mycotoxins using LC-MS/MS instrumentation.