SPE-LV6

Compact solid phase extraction system

Make your sample preparation process more reliable and secure

Innovative technology for all applications



PRESENTATION

Compact solid phase extraction system periods Make your analytical process more reliable and secure

The SPE-LV6 is a fully automated positive pressure workstation for sample preparation by SPE solid phase extraction. Specially developed for all analytical chemists looking for performance and alternatives to manual use, the SPE LV6 ensures method development and routine applications are executed with ease. Equipped with an innovative and robust technology, the SPE LV6 saves time and increases productivity thanks to the full automated process. The SPE LV6 uses calibrated tubing to grant high reproducibility and efficiency.

Certified LVD & EMC International Directives / CE / UL / CSA Including national differences for EU Group Differences, USA, Canada



Compatible with:

- SPE columns 1mL.3mL.6mL
- FPLC column 1mL

APPLICATIONS

- Environmental (waters, soils, pesticides...)
- Food & Feed (pesticides, drugs...)
- Biological / Toxicology (blood, urine, protein & peptide purification...)
- Petrochemical (Aromatic hydrocarbon...)





BENEFITS FOR YOUR SAMPLE PREPARATION

- 1. Fully automated system ensures time savings for users during sample preparation
- 2. Guarantees repeatable results
- 3. Create your purification methods with flexibility & simplicity
- 4. Improve your sample purification step
- 5. Automatic rinsing step avoids cross contamination
- 6. Innovative & robust technology
- 7. Plug and play in few minutes
- 8. Easy to use without specific knowledge
- 9. Cleaner & safer instrument with solvent waste recovery

SAFER & CLEANER

Keeps a cleaner and safer laboratory with solvent waste recovery Automatic solvent waste separations





Organic




Advion Interchim



7


TECHNOLOGY



POSITIVE PRESSURE TECHNOLOGY

Efficient & repeatable solvent & sample injections

- Constant flow rate from 50µL/min up to 15mL/min
- Increase molecule recovery
- Better separation & retention
- Avoids drying sorbent
- Allows work on resin & ion exchange sorbent








SAMPLE PURIFICATION METHOLOGY


SAMPLE VOLUME RANGE

Removes impurities Analyse molecule of interest



Load & clean your sample Isolate your molecule

ImpuritiesMolecule





For small volume

From 0,1mL Up to 25mL



Autosampler With sample loop





For large volume
From 10mL Up to 5L










REPRODUCIBILITY


Injection mode Sample volume		Autosampler (with needle) + sample loop			Large volume (syringe pump)	
		1 mL	3mL	5mL	50mL	
Sampling	RSD (%) Accuracy	0,33 0,991 ± 0,003mL	0,04 2,994 ± 0,001mL	0,04 4,983 ± 0,002mL	X	
Collection	RSD (%) Accuracy	0,75 0,96 ± 0,01mL	0,26 2,89 ± 0,01 mL	0,19 4,82 ± 0,01mL	0,14 49,67 ± 0,07mL	

ZERO CROSS CONTAMINATION

Thanks to automatic rinsing steps and dead-volume control

••










REDUCE CONCENTRATION/ EVAPORATION TIME

Drying step efficiency up to 92% in 30 seconds




ADAPTER

Universal adapters for SPE columns 1mL, 3mL & 6mL Creates a perfectly tight seal











INTEGRATED SOFTWARE

- Intuitive software
- 9" touch screen
- Create & Ssave hundred methods
- Start method in few seconds
- Automated report recording
- Method & data export with USB key
- Offline version to work on a laptop
- Start a run with a mixed column sizes & methods to optimize R&D





COMPACT FOOTPRINTS



Easy to set up

Installation on a bench or fume hood

No technical requirements for the installation

••

Advion Interchim



INSTRUMENT & ACCESSORIES -**PART NUMBERS**

Part Number	Description	Qty/Pack
SPE000	Interchim SPE-LV6 Instrument	1
SPE010	Sample Rack 13x100mm, 9mL	1
SPE020	Sample Rack 18x150mm, 25mL	1
SPE030	Collection Rack Vials 12x32mm, 1.8mL	1
SPEO40	Collection Rack Vials 12x75mm, 5mL	1
SPE050	Collection Rack Vials 13x100mm, 9mL	1
SPEO60	Collection Rack Vials 16x10mm, 11mL	1
SPEO70	Collection Rack Vials 18x150mm, 25mL	1
SPE300	Collection Rack Vials 28x140mm, 55mL	1



INSTRUMENT & ACCESSORIES -PART NUMBERS

Part Number	Description	Qty/Pack
SPE110	Insert Adapter SPE Columns, 3mL	1
SPE120	Insert Adapter SPE Columns, 1 mL	1
SPE510	Insert Adapter FPLC Columns, 1 mL	1
SPE520	SPE Columns 1 mL - 3 mL - 6 mL Adapter, PP	1
SPE530	SPE Columns 1 mL - 3 mL - 6 mL Adapter, PEEK	1
SPE540	Tubing 23cm Kit for Adapter & FPLC Column	1
468452	Universal PEEK Uptisur fitting 1/16"	1











CANISTERS & GLASSWARES - PART NUMBERS

Part Number	Description	Qty/Pack
CLEAR GLASS BO	ITILES	
B689C0	Bottle GL45, 100mL Clear Glass Duran Graduated with screw Cap and Pouring Ring PP	10
B689C1	Bottle GL45, 100mL Clear Glass Duran Graduated with screw Cap and Pouring Ring PP	1
B689E0	Bottle GL45, 150mL Clear Glass Duran Graduated with screw Cap and Pouring Ring PP	10
B689E1	Bottle GL45, 150mL Clear Glass Duran Graduated with screw Cap and Pouring Ring PP	1
BV2762	Bottle GL45, 250mL Clear Glass Duran Round	10
BV2761	Bottle GL45, 250mL Clear Glass Duran Round	1
HP6192	Bottle GL45, 500mL Clear Glass Duran Round	10
HP6191	Bottle GL45, 500mL Clear Glass Duran Round	1
B4PQA1	Bottle GL45, 3500mL Clear Glass Duran Round	1
HP6165	Bottle GL45, 1L Clear Glass Duran Round	1
HP6201	Bottle GL45, 2L Clear Glass Duran Round	1
B4PQA1	Bottle GL45, 3500mL Clear Glass Duran Round	1
JI6915	Bottle GL45, 5L Clear Glass Duran Round	1
JM5535	Bottle GL45, 10L Clear Glass Duran Round	1
B689V0	Bottle GL45, 15L Clear Glass Duran Graduated with screw Cap and Pouring Ring PP	1
B689X0	Bottle GL45, 20L Clear Glass Duran Graduated with screw Cap and Pouring Ring PP	1
CANISTER PE-HD		
	Interchim SPE-LV6, Canister PE-HD S60/61, 5L	1
	Interchim SPE-LV6, Canister PE-HD S60/61,12L	1
	Interchim SPE-LV6, Canister PE-HD S60/61, 20L	1
	Interchim SPE-LV6, Canister PE-HD S60/61, 30L	1



TEST TUBES - PART NUMBERS

Part Number	Description	Qty/Pack
TEST TUBES & CAPS		
AXETN5	Test Tubes, 12,00x75mm Borosilicate Glass EXP51, 6mL Without RIM - Round Botttom - 0,5-0,6mm	4 x 250
	Test Tubes, 13,00x100mm Borosilicate Glass EXP51, 10mL Without RIM - Round Botttom - 0,5-	
AXETO5	0,6mm	4 × 250
AXETP5	Test Tubes, 16,00x100mm Borosilicate Glass EXP51, 5mL Without RIM - Round Botttom - 0,5-0,6mm	4 × 250
AXETT5	Test Tubes, 18,00x150mm Borosilicate Glass EXP51, 28,5mL Without RIM - Round Botttom - 0.8mm	2 × 250
341651	10mL Screw Vials, 13mm lear Glass 13x100mm Rounded Botttom (For Sample Rack SPE010)	1 x 125
	13mm Screw Caps White Silicone/REDPTFE 1.5mm - Cross Slit (For Sample Rack SPE010	
BF2240	& Test Tubes 341652)	1 x 100
ARX3C0	Culture Tubes, 13,00x100mm Borosilicate Glass EXP51, 10mL 13-415 - Round Botttom - 0.9mm	4 x 250











How to request more information, a quotation or to place an order

Please see from the list below:

E-mail

Online

info.EU@advion-interchim.com quotes.EU@advion-interchim.com orders.EU@advion-interchim.com www.advion-interchim.com

Phone

+33 4 70 03 88 55

All your orders will be recorded & processed as quickly as possible.

You need technical assistance?
Our scientific experts are here to help.

instrumentation@advion-interchim.com

Advion Interchim Scientific is the market-facing name of a partnership between Advion, Inc. and Interchim SAS, separate companies that have agreed to market and distribute their products together.

All registered and trademark belong to their respective owners. Non-contractual pictures.