



Determination of Polycyclic Aromatic Hydrocarbons in Drinking Water by Liquid-Solid Extraction and HPLC with Coupled Ultraviolet and Fluorescence Detection*

UCT Part Numbers:

ECUNIPAH (2000 mg unendcapped C18, 83 mL cartridge)

Or

EUC1812M15 (2000 mg unendcapped C18, 15 mL cartridge)

CUC181M6 (1000 mg unendcapped C18, 6 mL cartridge)

ECSS25K (Anhydrous Sodium Sulfate)

EPA Method 550.1

Procedure

1. Cartridge Preparation

- a) Wash with 4 x 10 mL aliquots of methylene chloride (MeCl_2)
- b) Wash with 4 x 10 mL aliquots of methanol (MeOH)
- c) Wash with 2 x 10 mL aliquots of reagent water

Do not let the cartridge dry out after step 1) c. otherwise repeat starting at 1) b.

2. Sample Extraction

- a) Adjust the vacuum setting for a flow rate of 10-15 mL per minute
- b) Add the 1 liter sample to the cartridge
- c) Rinse sample bottle with reagent water, add to cartridge and draw through
- d) Dry cartridge by drawing full vacuum for 10 minutes

3. Sample Elution and Drying

- a) Elute the cartridge dropwise by using 2 x 5 mL aliquots of MeCl_2 and collect
- b) Rinse sample container with 5 mL of MeCl_2 , add to cartridge and draw through
- c) Prepare a drying column/funnel containing 10-20 g sodium sulfate by rinsing with 10 mL of MeCl_2 and discard
- d) Add the eluate to the drying column, draw through and collect
- e) Rinse the eluate vial and drying column with a 2 x 5 mL aliquot of MeCl_2 and collect

4. Sample Evaporation

- a) Evaporate the extract using a gentle stream of N₂ with a water bath or heating block temperature of 40°C. Evaporate to about 1.0 mL
- b) Add 3.0 mL of acetonitrile (ACN)
- c) Concentrate to a final volume of 0.5 mL

5. Sample Analysis

- a) Inject 5 - 100 µL into an HPLC

*See "Determination of Polycyclic Aromatic Hydrocarbons in Drinking Water by Liquid-Solid Extraction and HPLC with Coupled Ultraviolet and Fluorescence Detection", W. J. Bashe & T.V. Baker (Technology Applications, Inc, Environmental Monitoring Systems Laboratory, US Environmental Protection Agency, Cincinnati, OH

DCN-017009-191