



THC-COOH FROM ORAL FLUIDS BY LC-MS/MS OR GC-MS CLEAN SCREEN[®] DAU EXTRACTION COLUMN

Part #

ZSTHC020 - CLEAN SCREEN[®] THC 200 mg, 10 mL Tube

or

CSDAU206 - CLEAN SCREEN[®] THC 200 mg, 6 mL Tube

SMSTFA-1-1 – SELECTRA-SIL[®] MSTFA w/ 1% TMCS

SBSTFA-1-1 – SELECTRA-SIL[®] BSTFA w/ 1% TMCS

SLDA50ID21-5UM – Selectra[®] DA HPLC Column, 50 x 2.1 mm, 5 μm

1. PREPARE SAMPLE:

To 1 mL of oral fluid specimen 50ng/mL internal standard (THCA-D9) and let sit for ten minutes at room temperature

Mix/vortex for 10 seconds

Add 0.5 mL of Glacial Acetic Acid and vortex for 10 seconds

2. CONDITION CLEAN SCREEN[®] EXTRACTION COLUMN:

1 x 3 mL CH₃OH.

1 x 3 mL D.I. H₂O.

1 x 1 mL 0.1N HCl.

NOTE: Aspirate at full vacuum or pressure

3. APPLY SAMPLE:

Load at 1 to 2 mL/minute.

4. WASH COLUMN:

1 x 2 mL D.I. H₂O

1 x 2 mL 0.1 M HCl: Acetonitrile (70:30)

Dry column (10 minutes at full vacuum or pressure)

1 x 200 μL Hexane

Aspirate at full vacuum or pressure (Additional step to remove any residual moisture)

5. ELUTE ANALYTE:

1 x 2 mL Hexane/ Ethyl Acetate (50:50)

Collect eluate at 1 to 2 mL/minute

NOTE: Before proceeding, insure there are no water droplets at the bottom of the collection tube. This may increase drying time and decrease MSTFA derivatizing efficiency

6. DRY ELUATE:

Evaporate to dryness at < 40 °C.

7. RECONSTITUTE / DERIVATIZE:

- **LC-MS/MS:** Reconstitute sample in 100 μL of mobile phase
Inject 20 μL.
- **GC-MS:** Dissolve residue in 50 μL of Ethyl Acetate and 50 μL MSTFA w/
1% TMCS
Overlay with N₂ and cap. Mix/vortex
React 30 minutes at 70°C; Cool and inject 1 μL


Alternate Derivatization

1. Form TMS Derivatives by adding 50 μL BSTFA w/ 1% TMCS and 50 μL of Ethyl Acetate; React 45 minutes at 70 °C

Contributed by:

Janet Putnam, Assistant Laboratory Director/RP Advanced Toxicology Network, Memphis, TN

INSTRUMENT CONDITIONS (LC-MS/MS):



Analyte	MRM Transitions		Relative Retention Time (minutes)
	Q1	Q3	
1. THC-DELTA 9-COOH D ₉	352	308	2.44
2. THC-DELTA 9-COOH	343	299	2.49

PARAMETERS

Mobile Phase A: 0.1% Formic Acid in D.I. H₂O

Mobile Phase B: 0.1% Formic Acid in Methanol

Flow Rate: 0.5 mL/minute

Polarity: Negative

Reconstitute: 100 µL

Injection Volume: 20 µL


LC Column: Selectra[®] DA HPLC Column 50 x 2.1 mm 5 µm

Instrument: API 3200 Qtrap MS/MS with Shimadzu Prominence UFLC

Gradient:

Time	%A	%B
0.00	60	40
2.00	30	70
2.50	10	90
2.51	60	40
4.00	STOP	

INSTRUMENT CONDITIONS (GC-MS):



MSTFA/BSTFA TMS IONS

Analyte	Quantify Ion	Qualifier Ion 1	Qualifier Ion 2	Relative Retention Time (min)
1. THC-COOH	371	473	488	7.34
THC-COOH D ₃	374	476	491	7.31

PARAMETERS

GC/MS: Thermo ISQ Trace 1300

GC capillary column: 30 m x 0.25 mm (0.25 µm) TG-1MS

Injector: 1 µL Splitless, 250 °C

Oven temperature program: 170 °C (1) to 310 °C (30 °C/ minute): hold (5 minutes)

Carrier gas: Helium (1.2 mL/ minute)

MSD condition: Aux temperature: 280 °C, MS Source: 350 °C, MS Quad: 150 °C