



BATH SALTS IN BLOOD, PLASMA/SERUM, URINE, OR TISSUE BY LC-MS/MS OR GC-MS CLEAN SCREEN XCEL[®] I EXTRACTION COLUMN

Part #:

CSXCE111 – CLEAN SCREEN XCEL[®] 130 mg, 1 mL Tube

PFAA-0-1 – SELECTRA-SIL[®] PFAA

SPFPOH-1 – SELECTRA-SIL[®] PFPOH

SLDA50ID21-5UM – SELECTRA[®] DA HPLC Column, 50 x 2.1 mm, 5 μ m

1. PREPARE SAMPLE

To 1 mL of 100 mM phosphate buffer (pH 6.0) add internal standards
Add 1 -2 mL of blood, plasma/ serum, urine, or 1 g (1:4) tissue homogenate
Mix/vortex and let stand for 5 minutes
Add 2 mL of 100 mM phosphate buffer (pH 6.0). Mix/vortex
Sample pH should be 6.0 \pm 0.5.
Adjust pH accordingly with 100 mM monobasic or dibasic sodium phosphate.
Centrifuge for 10 minutes at 2000 rpm and discard pellet

2. APPLY SAMPLE

Load sample directly to column without any preconditioning.
Pull sample through at a rate of 1-2 mL/ minute.
Dry column thoroughly under full vacuum or positive pressure for 1 minute.

3. WASH

1 x 3 mL 98% Methanol: 2% Acetic Acid
Dry column thoroughly under full vacuum or positive pressure for a minimum of 5 minutes.

4. ELUTION

1 x 3 mL CH₂Cl₂/ IPA/ NH₄OH (78:20:2)
Collect eluate at 1 to 2 mL/minute.

NOTE: Prepare elution solvent daily.
Add IPA/ NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

5. DRY ELUTE

Add 50 μ L of 1% HCl in CH₃OH to each tube
Evaporate fraction to complete dryness under stream of dry air or nitrogen at ~ 35 °C.


NOTE: A 1% HCl in CH₃OH solution has been used to prevent volatilization by the formation of the hydrochloric salt of the drugs.

6. RECONSTITUTE / DERIVATIZE

- **LC-MS/MS:** Reconstitute sample in 100 μ L of mobile phase
Inject 5 μ L.
- **GC-MS:** Fluoroacylate with PFPA (PFAA)
Add 50 μ L PFPA. Over lay with N₂ and cap
*Improved derivatization by addition of PFPOH
React 20 minutes at 70 °C. Evaporate to dryness <40 °C
Reconstitute with 100 μ L Ethyl Acetate

NOTES: (It is important to dry the column thoroughly to achieve the highest recovery of all compounds. Any residual moisture will slow down the drying of the elution solvents prior to derivatization for GC/MS analysis, if being used. Also, any residual moisture could reduce the reactivity of the derivatization agent resulting in low GC/MS sensitivity.)

INSTRUMENT CONDITIONS (LC-MS/MS):



Analyte	MRM Transitions		Relative Retention Time (minutes)
	Q1	Q3	
1. Flephedrone	182.1	164.2	1.41
2. Methylone	208.1	160.1	1.66
3. Methadrone	194.1	161.1	1.96
4. Methedrone	178.1	145.1	2.34
5. Methethcathinone	192.2	174.0	2.98
6. MDPV	276.2	126.1	7.34
7. Pyravalerone	246.2	105.2	11.24

PARAMETERS

Mobile Phase A: 0.1% Formic Acid in D.I. H₂O

Mobile Phase B: 0.1% Formic Acid in Methanol

Flow Rate: 0.7 mL/minute

Polarity: Positive

Reconstitute: 100 µL

Injection Volume: 5 µL


LC Column: Selectra[®] DA HPLC Column 100 x 2.1 mm 5 µm

Instrument: API 4000 Qtrap MS/MS with Agilent 1200 Binary Pump SL

Isocratic:

Time	%A	%B
0.00	70	30
15.00	STOP	

INSTRUMENT CONDITIONS (GC-MS):



Fluoroacrylate with PFAA (PFAA) ions

Analyte	Quantify Ion	Qualifier Ion 1	Qualifier Ion 2	Relative Retention Time (min)
1.Flephedrone	123	204	160	7.32
2.Mephedrone	204	160	149	7.99
3.MCAT	218	174	91	8.33
4.Methylone	353	204	160	8.74
5.Methadrone	135	160	204	9.17
6.Ethylone	218	190	367	9.47
7.Butylone	218	160	367	9.61
8.Pyralvalerone	126	84	91	9.72
9.MDPV	126	96	84	10.62

PARAMETERS

GC/MS: Thermo ISQ Trace 1300

GC capillary column: 30 m x 0.25 mm (0.25 μ m) TG-1MS

Injector: 1 μ L Splitless, 250°C

Oven temperature program: 50 °C (1) to 310 °C (25 °C/ minute): hold (3.6 minute)

Carrier gas: Helium (1.2 mL/ minute)

MSD condition: Aux temperature: 280 °C, MS Source: 250 °C, MS Quad: 150 °C

Reference:

Comprehensive Forensic Toxicological Analysis of Designer Drugs; NIJ Grant

Author(s): Anthony P. DeCaprio, W. Lee Hearn, Madeleine J. Swortwood

Document No.: 244233

Date Received: December 2013